



## Changing the Number of Decimal Places in the EnviroQuant / Enhanced Summary Report

### MS EnviroQuant G1032C and G1701AA Enhanced or Environmental ChemStation

The file EPA.MAC controls the number of decimal places in the EnviroQuant summary report. This file can be edited to give the precision the user needs.

**Before making any changes to system macro files, BACKUP the file first. Suggested names for a backup copy of this file: EPA.FAC (for factory), or EPA.OLD.**

Edit \hpchem\msmacros\epa.mac, using Notepad or Wordpad. Search for the 4 occurrences of **FORMAT\$=** as shown below, and change the **###.##** ( or **#####.##**) from the end of each of these lines to the number of decimal places you need. For example, to change from 2 decimal places (the default ) to 3 decimal places, change **###.##** to **###.###**. If your calibration range requires more than 3 decimal places, you should maintain the same field width by reducing the number of places to the left of the decimal. For example, to change to 4 places to the right of the decimal, change **###.##** to **##.####**.

The following example shows the original macro lines and surrounding code. The shaded characters can be replaced by the characters below the shaded area to print 3 places to the right of the decimal. The small portion of code that must be edited is shaded. Try to maintain the original field width (6 or 7 characters wide) as shown in this example. Note that there are several occurrences of lines containing **FORMAT\$=** that are NOT involved in this change. This example is from G1701AA.

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material. The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose.

# HPCHEM\MSMACROS\EPA.MAC

```
!          format$="3#) 26% ###.## ##### #####  ###.## 5% 5%"
!eqnhance added % after area format for 'm' flag
          format$="3#) 26% ###.## 4% #####%  ###.## 5% 5%"
                                     ###.##
                                     ###.##
```

```
time = cal_pk_ret/60000-exp_rt
      if abs(time)<0.01
        time = 0
      endif
      Sprint rrt$ using "##.##",time
      Pn$=peak_name$ !eqnhance length restriction
      if tofile=1
!          print using #2,"  "+format$+ "/" ,result_num-
1,Peak_name$[1:26],cal_pk_ret/timediv,scan_num,pk_area,cal_pk_amt,amt_units$[1:5],rrt$
!eqnhance next line
```

```
-----
Pn$=peak_name$ !eqnhance length restriction
if mi_used = 0
  sprint qvalue$ using "3#",qvalue
else
  qvalue$="  "
endif
if tofile=1
  format$=" 3#) 26% ###.## 4% #####%  ###.## 5%1%  3%/"
                                     ###.##
                                     ###.##
  print using
#2,format$,result_num,pn$[1:26],cal_pk_ret/timediv,scannum$,pk_area,amod$,cal_pk_amt,amt_units$[1:5
],mod$,qvalue$
endif
if quant_printer=1
  format$="3#) 26% ###.## 4% #####%  ###.## 5%1%  3%"
                                     ###.##
                                     ###.##
  sprint lt$ using
format$,result_num,pn$[1:26],cal_pk_ret/timediv,scannum$,pk_area,amod$,cal_pk_amt,amt_units$[1:5],m
od$,qvalue$
  strprint lt$,indent,compoundrow
  compoundrow=compoundrow+1
  epacheckend
endif
-----
```

```

endif

!      format$="3#) 26% ###.## #### #####  ###.## 5%7%"
!eqnhance format
      format$="3#) 26% ###.## 4% #####%  ###.## 5% 7%"
      ###.##
      pn$=peak_name$ !eqnhance length restriction
      if tofile=1
!      print using #2," "+format$+"/",result_num-
1,peak_name$[1:26],cal_pk_ret/timediv,scan_num,pk_area,cal_pk_amt,amt_units$[1:5],mod$
!eqnhance
      print using #2,"
"+format$+"/",p_res,pn$[1:26],cal_pk_ret/timediv,scannum$,pk_area,amod$,cal_pk_amt,amt_units$[1:5],
mod$
      print using
format$,p_res,pn$[1:26],cal_pk_ret/timediv,scannum$,pk_area,amod$,cal_pk_amt,amt_units$[1:5],qsmo$
QaLogSurrogate      ! dwp 8/18/95
      if samt > 0
          print using #2," %/",surr$
      endif
endif
if quant_printer=1

```