

## 7500 ICP-MS Supplies

### Gas Line Supplies

#### Gas Line Supplies

Description	Specs	Usage	Unit	Part No.	Price
Inner sleeve for carrier/ auxiliary/blend gas			5/pk	5042-0922	\$28
Inner sleeve for plasma gas line			5/pk	5042-0923	28
Argon tubing assembly, 5 m Teflon tubing, one in-line filter and connector set				G1820-65023	230
Plasma and auxiliary gas line tubing, Tygon	3.17 mm id, 6.35 mm od, 35 cm (for plasma gas line). 2.38 mm id, 3.96 mm od, 35 cm (for auxiliary gas line), 1 ea.	Connects the plasma gas and auxiliary gas lines to the torch.		G1833-65428	21
Connector for end cap and gas line				5063-5263	38
Tubing for carrier/ makeup gas line, PTFE	3 mm id, 4 mm od, 350 mm		2/pk	G1833-65414	88

# 7500 ICP-MS Supplies

## Connectors



Connectors for Plasma and Auxiliary Gases  
P/N G1820-65027



Gas Line Connector  
P/N G1820-65052



Carrier Gas Connector  
P/N G1820-65214



Blend Gas Line Connector  
P/N G1820-65119

### Connectors

Description	Specs	Usage	Part No.	Price
Connectors, Teflon (2/pk)	6 mm and 4 mm	For plasma and auxiliary gases	G1820-65027	\$84
Carrier gas connector for Babington Nebulizer	L-shaped Teflon connector with a ferrule and O-ring	For connection with 4 mm diameter tubing	G1820-65214	166
Blend gas line connector, polypropylene		Connects 4 mm tubing with Babington and Cross Flow Nebulizers	G1820-65119	162
Carrier gas line connector, polypropylene	Includes ferrules and Viton O-ring	Connects 4 mm carrier gas tubing to Concentric Nebulizer	G1820-65052	104