



# Agilent Technologies

Innovating the HP Way

## Part number for GC/MS Interface and related parts Applies to 5973A/N MSD

### Part numbers for 5973A/N GC/MS Interface and related parts:

GC/MSD Interface ( <b>Complete</b> ) <b>(Includes heater/sensor and heater sleeve)</b>	G1099-60300
Welded Interface assembly <b>(Machined Interface only)</b>	G1099-60301
GC/MSD Interface O-Ring	0905-1405
GC/MSD Heater Sleeve <b>(Heater clamp)</b>	G1099-20210
Interface Heater/sensor assembly	G1099-60107
Interface Insulation	G1099-20301
Interface cover	G1099-00005
Interface/Cover mounting screws	0515-0380
Interface Column nut	05988-20066
Heater/Sensor set screw	0515-0236
Heater Sleeve screws <b>(Size M4 x 0.7)</b>	0515-0383

This document is believed to be accurate and up-to-date. However, Agilent Technologies, Inc. cannot assume responsibility for the use of this material.  
The information contained herein is intended for use by informed individuals who can and must determine its fitness for their purpose.