

G2350A Front assemblies




Figure 158. Front assemblies

Number	Description	Part number	Qty/Assy
1	Thumb nut	0535-0165	1
2	Gas drawer disconnect brazement	G2350-80680	1
3	O-ring	0905-1162	3
4	Emission source assembly	G2350-60560	1
5	Water drawer	G2350-61250	1
6	Spectrometer assembly	G2350-61170	1
7	Gas drawer	G2350-60530	1
NS	Ship kit	G2350-60800	-
NS	AED II Packaging assembly	G2350-61140	-