

4500 ICP-MS Supplies

Miscellaneous Components

Miscellaneous Components

Description	Specs	Usage/Replacement	Unit	Part No.	Price
O-ring for vacuum release	Viton O-ring	Used to seal between the vacuum chamber and its release valve	5/pk	G1820-65032	\$14
O-rings for vacuum chamber	Viton O-ring (one ring and one string type)	Used to seal between the vacuum chamber and its cover	2 ea.	G1820-65031	16
Filter for cooling water	One filter			G1820-65018	19
Filter kit for cooling water	Includes: filter for cooling water, housing for the water filter, connector for the water filter, stand for the water filter, and tubing for cooling water.			G1820-65017	570
Tool kit	Includes: nut driver (1/4 in., 10 mm); hexagonal wrench (H1.5, H3, H6); flat-blade screwdriver Phillips screw-driver set #1 and #2; short Phillips screwdriver, TV screwdriver, wrench for mounting the turbo with bolts; open end wrench set (H11x13, H10, H7/16, H9/16, H4, H2.5, and H5.5); tweezers, interface wrench, Babington Nebulizer removal tool, Babington Nebulizer sample line tool.			G1820-65019	408
Air filter		For side panel		G1820-65074	69
Air filter		For bottom		G1820-65073	69