

# 7500 ICP-MS Supplies

## Peristaltic Pump Tubing and Supplies

### *Peristaltic Pump Tubing and Supplies*

Description	Specs	Usage/Replacement	Unit	Part No.	Price
Tygon tubing	1.02 mm id, length = 400 mm	Standard tubing for sample introduction	12/pk	G1820-65215	\$99
Ismaprene tubing	1.52 mm id, length = 400 mm	Standard tubing for drainage; replacement required when tubing becomes stretched	12/pk	G1820-65216	99
ISTD tubing, Tygon	0.19 mm id, length = 400 mm	Used for on-line internal standard addition. Replacement is required when tubing is stretched	12/pk	G1820-65220	87
On-line ISTD addition Kit		Permits mixing of sample solution with internal standard solution at a fixed rate using peristaltic pump		G1820-65479	499
Drain tubing kit		Used to drain spray chamber to waste		G1820-65221	123
Drain tank				G1820-65016	215
Drain tank for organic solvent				G1820-65505	179
Tubing for drainage	4 mm id, 2.5 m			G1820-65515	96

#### **Helpful Hints:**

- Use Agilent's drain tubing kit to replace a heavily contaminated drain tubing assembly or when pumping efficiency is low.
- Leaving the tubing attached to the pump will shorten the lifetime of the tubing. It is recommended to release the peristaltic pump tubing when the system is not in use.